

## INFORMATION BULLETIN



# THE ROLE OF PREANALYTICAL FACTORS IN CHEMISTRY AND IMMUNOASSAY TESTING.

## INTRODUCTION >

Measurement of biochemical markers is an important aid to clinicians in the early detection, diagnosis, monitoring and prognosis of disease. Specimen quality plays a key role in ensuring accuracy of those measurements in clinical laboratory testing. To gain efficiencies in workflow and decrease turnaround time (TAT), many laboratories have adopted new strategies and practices, including transitioning from:

- › glass to plastic specimen collection tubes
- › serum to anticoagulated plasma samples
- › manual processing to lab automation
- › sample collection by laboratory staff to non-laboratory personnel

As laboratories automate more processes, less time is dedicated to sample inspection steps and, consequently, monitoring specimen quality. Preanalytical factors can be magnified by sensitive immunoassays and present an increasing challenge to quality clinical care.

## PREANALYTICAL VARIABLES THAT COULD AFFECT RESULTS >

As much as 84 percent of laboratory errors can be attributed to the preanalytical phase of clinical laboratory testing. This phase is comprised of an evaluation of patient condition, as well as specimen collection, transport, processing and placement on the analyzer.<sup>1,2,3,4,5</sup> Patient samples with circulating protein interferants, such as human anti-animal antibodies (HAAA) and rheumatoid factor (RF), may affect the results of certain assays. These are also examples of potential sources of error outside the control of the laboratory.<sup>6,7</sup> Knowledge of such factors is important when determining the appropriate interpretations of results.

The large majority of preanalytical errors are due to compromised sample quality<sup>1,2,3,4,5</sup> which can occur during specimen collection, storage, transport and processing. Common factors contributing to error include: incorrect labeling of tubes, insufficient blood draw volume, insufficient mixing, cellular contamination in plasma specimens and inadequate clotting of serum specimens.

To maintain sample quality, each stage in sample preparation is important. It is critical that personnel performing blood collection adhere to all recommendations specified by blood collection tube manufacturers. Deviations from the manufacturer's recommendations must be validated in individual laboratories.

## FACTORS AFFECTING PLASMA SAMPLES >

---

While serum may provide the cleanest sample from an interference perspective, certain factors can affect processing the sample in a timely manner. Because urgent, critical decisions are based on STAT results, heparinized plasma samples have become the preferred sample type and are widely used. Laboratory Medicine Practice Guidelines, published by the National Academy of Clinical Biochemistry (NACB), recommend plasma for STAT analysis of cardiac markers.<sup>8</sup> Plasma provides the best opportunity for achieving desired rapid turnaround time. However, there are variables that must be controlled to obtain the best possible sample for analysis.

Given that plasma samples contain anticoagulants, the cellular components (i.e., white blood cells, red blood cells and platelets) are not trapped in a clot during the normal coagulation process of a serum sample. Following centrifugation, plasma samples may still contain trace amounts of cellular material, as well as latent fibrin. Gel separator tubes reduce the incidence of resuspension of these formed elements. However, some materials (e.g. platelets) may remain above the plasma gel interface barrier. These factors can cause non-specific binding to the solid phase - the microparticles - leading to erroneous results.

## HEPARIN AS AN ANTICOAGULANT >

---

Heparin, a negatively-charged molecule used to inhibit clotting, can bind to some analytes, antibodies and cellular material. This can interfere with the antigen-antibody interaction in the test method.<sup>9,10</sup>

If a tube has insufficient blood volume, there is an excess of heparin. Maintaining an optimum, sample-to-additive ratio is important for effective anticoagulation and accurate laboratory tests.<sup>11,12,13</sup> A key step in the sample handling process is ensuring that the blood draw sample volume is at least 90 percent of the stated volume on the collection tube.<sup>11</sup> The Clinical and Laboratory Standards Institute (CLSI) has published guidelines for blood specimen handling. Heparin is also a commonly used pharmaceutical agent to inhibit clotting in critical care patients. Inadequate cleaning of an intravenous line prior to blood collection can also create an excess in the sample.

## POSSIBLE MECHANISMS THAT COULD INTERFERE WITH HEPARIN ANTICOAGULANT ACTIVITY

Certain mechanisms might interfere with heparin anticoagulant action, resulting in fibrin formation in a plasma sample.<sup>12</sup> These include:

1. **The ability of heparin to bind cell membranes/proteins, such as platelets.** Heparin has a tendency to bind to plasma proteins and cell membranes, thus making its pharmacological action unpredictable. The presence of cellular proteins and membranes could result in binding of heparin, therefore competing and interfering with anticoagulation.
2. **Upon re-exposure to heparin some patients will exhibit heparin-induced thrombocytopenia (HIT).** This condition can cause heparin-induced or heparin-facilitated platelet aggregation and resulting in platelet counts. The activated platelets release platelet factor 4 (PF4) that allows clotting by neutralizing heparin.

## EFFECT OF FIBRIN IN PLASMA AND SERUM SAMPLES >

---

Chemistry and immunoassays are susceptible to interference by fibrin. Small amounts of fibrin (and other membrane fragments or cell stroma) may affect sensitive immunoassays and chemistry tests. The presence of gross amounts of fibrin in the specimen (serum or plasma) may cause blockage of instrument sample aspiration probes. Such a blockage could lead to erroneous assay results.

### PLASMA SAMPLES

Inadequate collection-tube mixing may result in uneven distribution of the heparin additive throughout the specimen. This could lead to localized areas within the specimen where the anti-thrombin effect of the heparin is insufficient to prevent the formation of fibrin. Thus, it is essential to ensure thorough mixing by gentle inversion immediately after blood is drawn in the tube (following manufacturer's guidelines for each tube type). In the past a liquid anticoagulant was used in many glass tubes, facilitating easy mixing. Today the walls of the tube are coated with a powdered anticoagulant. This is not as easily mixed in the sample unless the required mixing occurs immediately after collection.

Since the heparin additive in specimens typically degrades over time, residual thrombin in the specimen can convert soluble fibrinogen to insoluble fibrin. Flocculent matter can frequently be observed in stored samples. Care should be taken to recentrifuge such samples prior to analysis.

### SERUM SAMPLES

Inadequate clotting time, improper mixing, and failure to place the tube in an upright position can lead to incomplete clot formation. Following centrifugation, the resulting sample may appear satisfactory with a defined layer of cells at the base of the tube and a clear layer of serum above. Despite this appearance, the clotting process may not have been completed prior to transportation, centrifugation and placement of the specimen on the analyzer. Further coagulation in the serum may subsequently occur, leading to the production of "latent" fibrin, which can interfere with the quality of a result.

For plastic tubes, thorough mixing by gentle inversion is essential to ensure even distribution of the clot activator throughout the specimen. This will also allow completion of the clotting process. Some cardiac patients will have therapeutic levels of anticoagulant in their blood. This could increase clotting time in the tube and thus increase the potential for the formation of "latent" fibrin in the preanalytical phase.

### CONCLUSION

Considering all of the above factors, serum appears to be the superior sample for immunoassays. Many laboratories use heparinized plasma for faster test turnaround times, and to avoid prolonged clotting times in patients with high circulating levels of heparin. Regardless of which sample type is used, following the blood collection tube manufacturer's specimen collection and handling recommendations will help to reduce preanalytical laboratory error. In order to minimize laboratory error due to specimen quality, the key preanalytical actions are:

1. Adequately fill the collection tube to the full volume
2. Ensure proper mixing immediately after collection
3. Allow adequate clotting time (minimum of 30 minutes) for serum specimens
4. Ensure proper centrifugation
5. Avoid resuspension of separated samples, including tubes with a gel barrier

## REFERENCES >

---

1. Plebani M, Carraro P. Mistakes in a stat laboratory: types and frequency. *Clin Chem* 1997;43:1348-51.
2. Stahl M, Lund ED, Brandslund I. Reasons for a laboratory's inability to report results for requested analytical tests. *Clin Chem* 1998;44:2195-97.
3. Wiwankit V. Types and frequency of preanalytical mistakes in the first Thai ISO 9002: 1994 certified clinical laboratory, a 6-month monitoring. *BMC Clin Path* 2001;1:5-9.
4. Bonini P, Plebani M, Ceriotti F, Rubboli, F. Errors in laboratory medicine. *Clin Chem* 2002;48:691-98.
5. Kaltra J. Medical Errors: impact on clinical laboratories and other critical areas. *Clin Biochem* 2004;37:1052-62.
6. Marks V. False-positive immunoassays results; a multicenter survey of erroneous immunoassay results from assays of 74 analytes in 10 donors from 66 laboratories in seven countries. *Clin Chem* 2002;48:2008-16.
7. Ismail et al. Wrong biochemistry results: two case reports and observational study in 5310 patients on potentially misleading thyroid-stimulating hormone and gonadotropin immunoassay results. *Clin Chem* 2002;48:2023-29.
8. Christenson R.H. et al. NACB Biomarkers of Acute Coronary Syndrome and Heart Failure (draft guidelines), 2004.
9. Zaninotto M, et al. Quality specifications for biochemical markers of myocardial injury. *Clinical Chimica Acta* 2004;346:65-74.
10. Gerhardt W, et al. Troponin T and I assays show decreased concentrations in heparin plasma compared with serum: lower recoveries in early than in late phases of myocardial injury. *Clin Chem* 2000;46:817-821.
11. Clinical and Laboratory Standards Institute. H1-A6; Tubes and Additives for Venous Blood and Capillary Blood Specimen Collection; Approved Standard-Sixth Edition 2010.
12. Bush V. Why doesn't my heparinized plasma specimen remain anticoagulated? *BD Lab Notes*. Spring 2003;13(2):9-12.
13. Dubrowny N. Smith S. Hanna S. From vein to the analyzer: sample handling impact. *BD-BCI workshop*. 2005:AACC annual meeting.



Beckman Coulter and the stylized logo are trademarks of Beckman Coulter, Inc. and are registered with the USPTO.

**For Beckman Coulter's worldwide office locations and phone numbers, please visit [www.beckmancoulter.com/contact](http://www.beckmancoulter.com/contact)**