


IMMUNOBLOT

Test Principle

Antigens are transferred to a nitrocellulose membrane using a micro-dispensing method.



User Comfort

- Ready-to-use components
- Colour-coded strips
- Interchangeable components
- Positive and Negative controls
- Control line on the strip
- Easy assay procedure

Advantages

- Identical assay procedure
- Easy interpretation and reproducibility of results
- High diagnostic specificity and sensitivity
- Ready for automation
- Customer support

Clinical Application

BLOT-LINE ANA

- Confirmatory test for EIA ENA screen plus
- Differential diagnosis of systemic autoimmune diseases by determination of specific ANAs

Protocol Summary

| Step | Test steps |
|---|---|
|  1. | Pipette Universal solution 2 ml |
|  2. | Strips soaking 10 min. at room temperature - Shaker |
|  3. | Aspirate |
|  4. | Dilute samples - serum/plasma 1:51 (30 µl + 1,5 ml) |
|  5. | Pipette Controls and diluted samples 1.5 ml |
|  6. | Incubate 30 min. at room temperature - Shaker |
|  7. | Aspirate samples and wash strips with 1.5 ml of Universal solution 3-times for 5 min. - Shaker |
|  8. | Pipette Conjugate 1.5 ml |
|  9. | Incubate 30 min. at room temperature - Shaker |
|  10. | Aspirate Conjugate and wash strips with 1.5 ml of Universal solution 3-times for 5 min. - Shaker |
|  11. | Pipette Substrate solution (BCIP/NBT) 1.5 ml |
|  12. | Incubate 15 min. at room temperature - Shaker |
|  13. | Aspirate Substrate solution and wash strips with 2 ml of distilled water 2-times for 5 min. - Shaker |
|  14. | Sticking and evaluation of strips |

Antigens


BLOT-LINE ANA

Recombinant antigens:

Ro52/SS-A, Ro60/SS-A, La/SS-B, RNP-A, RNP-C, RNP 68, SmB, SmD, Scl 70, Jo-1, Centromere B, Centromere A, PMScl, PL-7, PL-12, ribosomal protein P0, PCNA, dsDNA, Histones, M2, DFS70 and Nucleosome

Test Characteristics

| Immunoblot | Diagnostic Sensitivity | Diagnostic Specificity |
|---------------|------------------------|------------------------|
| BLOT-LINE ANA | 97.2% | 99.1% |



Clinical Data

Detection of individual antigens for group of patients with SLE - results of BLOT-LINE ANA

