

Classic Coagulation state of the art



Classic fully automated coagulation analyser

- New Design
- Increased Throughput
- Full Traceability









The success story of the Compact Series started in 1993. Now, more than 2000 systems are still being used throughout the world. A large part of its success is due to its compact form and reliability. In 2009 we have completely revamped the Compact Series.

The Compact is the classic coagulation analyser. Simple, Fast, Reliable.

The Compact X expands the measuring system to include Chromogenic and Immunologic tests such as AT III and D-Dimer. In addition, the intelligent cuvette bar return system ensures all cuvettes are used to their full capacity.

The advantages of the new generation:

- Increased throughput
- Open working area for more effective handling
- New software with comprehensive traceability – including reaction curves, reagent level display, QC database and the usual, intuitive user-friendly interface.

By selecting a Compact Series instrument, you are choosing a next-generation, fully automated coagulation analyser which combines all the advantages of Behnk Thrombolyzer reliability and modern analysis techniques.

to their full capacity.	sis techniques.	
Specification	Thrombolyzer Compact	Thrombolyzer Compact X
Sample throughput	160 tests/h	160 tests/h
Open system	for use with almost all reagents	for use with almost all reagents
On-board sample capacity	32	32
Reagent positions (cooled)	16	16
Cuvettes on-board	240 tests	240 tests
Reloadable at any time: Cuvettes, Samples, Reagents	√	✓
Automatic test repetition and reflex tests	✓	✓
Automatic calibration curve creation	√	√
Automatic cuvette rack return	-	√
Chromogenic and Immunologic tests	-	√
Measuring channels	4 at 620 nm	4 at 405 nm & 620 nm
Operating system	LINUX	LINUX
LIS Interface	RS 232 or LAN	RS 232 or LAN
Barcode scanner	sold seperately	sold seperately
Dimensions (LxWxH)	56 x 55 x 38 cm	56 x 55 x 38 cm
Weight	27,5 kg	27,5 kg
Sound Level	48 – 61 dB	48 – 61 dB
Mains voltage	100 – 240 V	100 – 240 V
Frequency range	47 – 63 Hz	47 – 63 Hz
Power consumption	100 VA / 0.34 BTU/h	100 VA/0.34 BTU/h
Supply range Thrombolyzer	Compact/Compact X, instrument ac	ccessories, computer, flat screen monito



Behnk Elektronik

Behnk Elektronik GmbH & Co. KG

Hans-Böckler-Ring 27 22851 Norderstedt

Germany

Telephone +49 (0) 40 - 529 86 10 Telefax +49 (0) 40 - 524 10 94

Web www.behnk.de E-mail info@behnk.de

Thrombolyzer Compact

Proven. Reliable. Fast.